

K2 Li CAD Glass Ceramic Bloc

Before and after the crystallization firing & customization



K2 Li CAD – High-strength Lithium Disilicate Glass Ceramic

K2 Li CAD blocks from Yeti Dental are the ideal solution for dental laboratories that refuse to compromise on aesthetics and efficiency. This modern glass ceramic was specially developed for digital workflows and combines natural light dynamics with exceptional stability.

Maximum Durability: High flexural strength (400 MPa) for long-lasting restorations
Natural Aesthetics: Optimized translucency for a lifelike play of light
Precise Processing: Excellent milling behavior without material chipping
Versatile Shade Range: Wide selection for precise color matching.
Open System: Fully compatible with common CAD/CAM systems
Dimensions: 18/15/13 mm

Indications: Single crowns (anterior and posterior), Veneers, Inlays and Onlays



Color	HT (High Translucency)	LT (Low Translucency)
A1	389-1011	389-1012
A2	389-1021	389-1022
A3	389-1031	389-1032
A3,5	389-1041	389-1042
B1	389-1061	389-1062
B2	389-1071	389-1072
C1	389-1101	389-1102
C2	389-1111	389-1112
D2	389-1141	389-1142
BL1	389-1171	389-1172
BL2	389-1181	389-1182
BL3	389-1191	389-1192
BL4	389-1201	389-1202