



Fixing Gel

Adhesive gel for secure fixation of ceramic applications on the model die.



Advantages of Fixing Gel:

- safe hold of the ceramic work on the model
- free layering
and check the articulation contacts,
without that the work is falling down on the floor
- thin bonding layer with very good bonding effect
- color-neutral and does not influence the ceramic color
- **burns without residue** -
no residues on the ceramic material
- water soluble -
with warm water or steam jet
easy to remove from the model

For which materials can Fixing Gel be used?

All-ceramics, glass-ceramics, zircon, zirconium oxide (ZrO₂), lithium disilicates and pressed ceramics, metal ceramics, non-precious alloys, precious metal alloys and dental alloys.

Fields of application:

Crowns, bridges, veneers, inlays.

